



Concise Separations取扱製品

株式会社エーディーエステック

Part #	Part Name	Body	Int. Dia.	Length
282013	Column Coupler, PEEK (connects columns to guard holders)	PEEK		
AXC-99-1300	Universal Holder For All 4.0x20mm Guard Cartridges	SS	4.0mm	24mm
<b>Amino Acid Analysis Guard Cartridges:</b>				
AAA-99-1311	AminoSep Lithium Cartridge 2/pk	Pktit	4.0mm	20mm
AAA-99-1312	AminoSep BG Sodium Cartridge 2/pk	Pktit	4.0mm	20mm
AAA-99-1353	AminoSep AA 911 Guard Cartridge 2/pk	Pktit	4.0mm	20mm
AAA-99-1354	AminoSep AA 511 Guard Cartridge 2/pk	Pktit	4.0mm	20mm
<b>Amino Acid Analysis Guard Kits (Includes Stainless Steel Holder + 2 Cartridges):</b>				
AAA-99-2311	AminoSep Li+ Guard Kit	SS/pktit	4.0mm	24mm
AAA-99-2312	AminoSep BG Na+ Guard Kit	SS/pktit	4.0mm	24mm
AAA-99-2353	AminoSep AA 911 Guard Kit	SS/pktit	4.0mm	20mm
AAA-99-2354	AminoSep AA 511 Guard Kit	SS/pktit	4.0mm	20mm
<b>Amino Acid Analysis Columns:</b>				
AAA-99-6310	AminoSep Beckman Gold Na+	SS	4.0mm	200mm
AAA-99-6311	AminoSep Beckman Style Li+	SS	4.0mm	100mm
AAA-99-6554	AminoSep Sodium AA 511HS	SS	4.6mm	120mm
AAA-99-7554	AminoSep Sodium AA 511	SS	4.6mm	150mm
AAA-99-8553	AminoSep Sodium AA 911	SS	4.6mm	250mm
AAA-99-8553W	AminoSep Sodium AA 911 Waters Compatible	SS	4.6mm	250mm
<b>Anion (Ion Chromatography) Analysis Guards:</b>				
ANX-99-0010	ICSep AN Guard Cartridge 3/pk (AN1 and AN300 compatible)	Peek	3.0mm	10mm
ANX-99-0014	ICSep AN1SC Guard Cartridge 3/pk	Peek	3.0mm	10mm
ANX-99-0015	ICSep AN2 Guard Cartridge 3/pk	Peek	3.0mm	10mm
ANX-99-0016	ICSep AN300B Guard Cartridge 3/pk	Peek	3.0mm	10mm
ANX-99-0018	ICSep AN2 HC Guard Cartridge 3/pk	Peek	3.0mm	10mm
ANX-99-0040	ICSep Holder for 3x10mm GC	Peek	3.0mm	10mm
ANX-99-3510	ICSep AN Guard Column (AN1 and AN300 compatible)	Peek	4.6mm	50mm
ANX-99-3514	ICSep AN1SC Guard Column	Peek	4.6mm	50mm
ANX-99-3515	ICSep AN2 Guard Column	Peek	4.6mm	50mm
ANX-99-3516	ICSep AN300B Guard Column	Peek	4.6mm	50mm
ANX-99-3517	ICSep AN2 HC Guard Column	Peek	4.6mm	50mm
AXC-99-0010	ICSep Column Replacement Frits 5/pk, PEEK collar	Peek	4.6mm	
AXC-99-0012	ICSep Column Replacement Frits 5/pk, Kel-F collar	Peek	4.6mm	
GRD-99-0703	ICSep AN Guard Disc Replacement, 5/pk	Membrane	4.6mm	
<b>Anion (Ion Chromatography) Analysis Columns:</b>				
ANX-99-3518	ICSep AN1HS	Peek	4.6mm	50mm
ANX-99-3566	ICSep AN300BHS	Peek	4.6mm	50mm
ANX-99-4650	ICSep AN SDA	Peek	4.6mm	50mm
ANX-99-5511	ICSep AN1 Fast Analysis	PEEK	4.6mm	100mm
ANX-99-5518S	ICSep AN1 SS	SS	4.6mm	100mm
ANX-99-7613	ICSep AN300H (for EPA Method 300.a)	Peek	5.5mm	150mm
ANX-99-7620	ICSep AN300 (for EPA Method 300.a)	Peek	5.5mm	150mm
ANX-99-7633	ICSep AN300 Column Kit (includes column, guard column, coupler)	Peek	5.5mm	150mm
ANX-99-7643	ICSep AN300 Cartridge Kit (includes column, guard holder, cartridges, coupler)	Peek	5.5mm	150mm
ANX-99-8505	ICSep AN2 Fast Analysis	Peek	4.6mm	150mm
ANX-99-8511	ICSep AN1	Peek	4.6mm	250mm
ANX-99-8514	ICSep AN1SC	Peek	4.6mm	250mm
ANX-99-8515	ICSep AN2	Peek	4.6mm	250mm
ANX-99-8516	ICSep AN300B (for EPA Method 300.b)	Peek	4.6mm	250mm
ANX-99-8518	ICSep AN2 HC, High Capacity	Peek	4.6mm	200mm
ANX-99-8531K	ICSep AN1 Column Kit (includes column, guard column, coupler)	Peek	4.6mm	250mm

ANX-99-8534K	ICSep AN1SC Column Kit (includes column, guard column, coupler)	Peek	4.6mm	250mm
ANX-99-8535K	ICSep AN2 Column Kit (includes column, guard column, coupler)	Peek	4.6mm	250mm
ANX-99-8536K	ICSep AN300B Column Kit (includes column, guard column, coupler)	Peek	4.6mm	250mm
ANX-99-8538K	ICSep AN2 HC Column Kit ( Includes column, guard column, coupler)	Peek	4.6mm	200mm
ANX-99-8541K	ICSep AN1 Cartridge Kit (includes column, guard holder, cartridges, coupler)	Peek	4.6mm	250mm
ANX-99-8544K	ICSep AN1SC Cartridge Kit (includes column, guard holder, cartridges, coupler)	Peek	4.6mm	250mm
ANX-99-8545K	ICSep AN2 Cartridge Kit (includes column, guard holder, cartridges, coupler)	Peek	4.6mm	250mm
ANX-99-8546K	ICSep AN300B Cartridge Kit (includes column, guard holder, cartridges, coupler)	Peek	4.6mm	250mm
ANX-99-8548K	ICSep AN2 HC Cartridge Kit ( includes column, guard cartridges, coupler, holder)	Peek	4.6mm	200mm
<b>Carbohydrate Analysis Guard Cartridges:</b>				
AXC-99-1300	Universal Holder For All 4.0x20mm Guard Cartridges	SS	4.0mm	24mm
CHO-99-1345	CarboSep CHO 413 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1351	CarboSep CHO 611 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1352	CarboSep CHO 611OH Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1353	CarboSep CHO 620 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1354	CarboSep CHO 682 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1355	CarboSep CHO 820 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1360	CarboSep CHO 87C Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1362	CarboSep CHO 87K Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1363	CarboSep CHO 87N Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1364	CarboSep CHO 87P Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1365	CarboSep CHO 87MM Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1367	CarboSep USP L19 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1371	CarboSep CHO 411 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1372	CarboSep CHO 782 Cartridge 2/pk	Pktit	4.0mm	20mm
CHO-99-1373	CarboSep CHO 882 Cartridge 2/pk	Pktit	4.0mm	20mm
<b>Carbohydrate Analysis Guard Kits Include Stainless Steel Holder + 2 Cartridges:</b>				
CHO-99-2345	CarboSep CHO 413 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2351	CarboSep CHO 611 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2352	CarboSep CHO 611OH Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2353	CarboSep CHO 620 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2354	CarboSep CHO 682 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2355	CarboSep CHO 820 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2360	CarboSep CHO 87C Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2362	CarboSep CHO 87K Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2363	CarboSep CHO 87N Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2364	CarboSep CHO 87P Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2365	CarboSep CHO 87MM Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2367	CarboSep USP L19 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2371	CarboSep CHO 411 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2372	CarboSep CHO 782 Guard Kit	SS/pktit	4.0mm	24mm
CHO-99-2373	CarboSep CHO 882 Guard Kit	SS/pktit	4.0mm	24mm
<b>Carbohydrate Analysis Guard Columns:</b>				
CHO-99-3560	CarboSep CHO 87C Guard Column	SS	4.6mm	50mm
CHO-99-3564	CarboSep CHO 87P Guard Column	SS	4.6mm	50mm
<b>Carbohydrate Analysis Columns:</b>				
CHO-99-5854	CarboSep CHO 682 USP 100, Lead Form	SS	7.8mm	100mm
CHO-99-5860	CarboSep CHO 87C Fast Analysis, Calcium Form	SS	7.8mm	100mm
CHO-99-5882	CarboSep CHO 882 Fast Analysis, Lead Form	SS	7.8mm	150mm
CHO-99-7752	CarboSep CHO 611OH , Sodium Form	SS	6.5mm	150mm
CHO-99-7770	CarboSep CHO 782, Lead Form	SS	7.8mm	300mm
CHO-99-8453	CarboSep USP L19, Calcium Form Mannitol/Sorbitol Analysis	SS	4mm	250mm
CHO-99-8770	CarboSep CHO 882, Lead Form	SS	7.8mm	300mm
CHO-99-8854	CarboSep CHO 682 200, Lead Form	SS	7.8mm	200mm
CHO-99-8855	CarboSep CHO 820, Calcium Form	SS	7.8mm	200mm
CHO-99-8864	CarboSep CHO 87P 200, Lead Form	SS	7.8mm	200mm
CHO-99-8865	CarboSep CHO 87MM, Mixed Mode Form	SS	7.8mm	200mm

CHO-99-9751	CarboSep CHO 611, Sodium Form Corn Syrup	SS	6.5mm	300mm
CHO-99-9753	CarboSep CHO 620, Calcium Form	SS	6.5mm	300mm
CHO-99-9753W	CarboSep CHO 620, Calcium Form Waters Compatible	SS	6.5mm	300mm
CHO-99-9850	CarboSep CHO 411, Sodium Form Oligosaccharides	SS	7.8mm	300mm
CHO-99-9853	CarboSep CHO 682 Fast Analysis, Lead Form	SS	7.8mm	150mm
CHO-99-9854	CarboSep CHO 682, Lead Form	SS	7.8mm	300mm
CHO-99-9855	CarboSep CHO 820, Calcium Form	SS	7.8mm	300mm
CHO-99-9860	CarboSep CHO 87C, Calcium Form	SS	7.8mm	300mm
CHO-99-9862	CarboSep CHO 87K, Potassium Form	SS	7.8mm	300mm
CHO-99-9863	CarboSep CHO 87N, Sodium Form	SS	7.8mm	300mm
CHO-99-9864	CarboSep CHO 87P, Lead Form	SS	7.8mm	300mm
CHO-99-9865	CarboSep CHO 87MM, Mixed Mode Form	SS	7.8mm	300mm
CHO-99-9866	CarboSep CHO 87P 100, Lead Form	SS	7.8mm	100mm
CHO-99-9867	CarboSep CHO 87P Fast Analysis, Lead Form	SS	7.8mm	150mm
<b>Carbohydrate Semi Prep Columns:</b>				
CHO-99-6353	CarboSep SemiPrep CHO 620, Calcium Form	SS	21.2mm	300mm
CHO-99-6360	CarboSep SemiPrep CHO 87C, Calcium Form	SS	21.2mm	300mm
CHO-99-6363	CarboSep SemiPrep CHO 87N, Sodium Form	SS	21.2mm	300mm
CHO-99-6370	CarboSep SemiPrep CHO 882, Lead Form	SS	21.2mm	300mm
<b>Cation Analysis:</b>				
CTX-99-3540	ICSep CN3 Metal Analysis Column	SS	4.6mm	40mm
CTX-99-1340	ICSep CN3 Guard Cartridges 2/pk	SS/pktit	4.0mm	24mm
CTX-99-2040	ICSep CN3 Guard Kit	SS/pktit	4.0mm	24mm
<b>Organic Acids Analysis Guard Kits (Includes Stainless Steel Holder + 2 Cartridges):</b>				
ICE-99-1372	Coregel Carbonate Analysis Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2352	Coregel USP L17 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2353	Coregel ARH 601 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2354	Coregel ORH 801 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2355	Coregel 87H Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2356	Coregel 107H Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2357	Coregel ION 310 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2358	Coregel ION 300 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2359	Coregel 801 FA Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2360	Coregel 64H Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-2361	Coregel 87H3 Guard Kit	SS/pktit	4.0mm	24mm
ICE-99-3510	Coregel WA1 Wine Analysis Guard Kit	SS/pktit	4.0mm	24mm
<b>Organic Acids Analysis Guard Column:</b>				
CHO-99-3561	Coregel 87H Guard Column (for 87H & 87H3)	SS	4.6mm	50mm
<b>Organic Acids Analysis Guard Cartridges:</b>				
AXC-99-1300	Universal Holder For All 4.0x20mm Guard Cartridges	SS	4.0mm	24mm
ICE-99-1310	Coregel WA1 Wine Analysis Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2362	Coregel USP L17 Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2363	Coregel ARH 601 Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2364	Coregel ORH 801 Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2366	Coregel 107H Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2367	Coregel ION 310 Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2368	Coregel ION 300 Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2369	Coregel 801 FA Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2370	Coregel 64H Cartridge 2/pk	Pktit	4.0mm	20mm
ICE-99-2371	Coregel 87H Cartridge 2/pk (for 87H & 87H3)	Pktit	4.0mm	20mm
ICE-99-2372	Coregel Carbonate Analysis Cartridge 2/pk	Pktit	4.0mm	20mm
<b>Organic Acids Analysis Columns:</b>				
ICE-99-5753	Coregel ARH 601 Aromatic Acids Analysis	SS	6.5mm	100mm
ICE-99-5861	Coregel 87H Fast Organic Acids Analysis	SS	7.8mm	100mm
ICE-99-7752	Coregel ION 310 Fast Organic Acids Analysis	SS	6.5mm	150mm

ICE-99-8461	Coregel USP L17	SS	4.0mm	250mm
ICE-99-9753	Coregel 801 Fast Analysis	SS	7.8mm	200mm
ICE-99-9754	Coregel ORH 801	SS	6.5mm	300mm
ICE-99-9810	Coregel WA1 Wine Analysis Column	SS	7.8mm	300mm
ICE-99-9850	Coregel ION 300	SS	7.8mm	300mm
ICE-99-9860	Coregel 64H	SS	7.8mm	300mm
ICE-99-9861	Coregel 87H3	SS	7.8mm	300mm
ICE-99-9865	Coregel 87H3 Fast Analysis	SS	7.8mm	150mm
ICE-99-9866	Coregel 107H	SS	7.8mm	300mm
ICE-99-9867	Coregel 107H Fast Analysis	SS	6.5mm	150mm
ICE-99-9871	Coregel 87H3, 6.5mm id	SS	6.5mm	300mm
ICE-99-9872	Coregel CA1 Carbonate Analysis	SS	7.8mm	250mm
<b>Organic Acid Semi Prep Columns:</b>				
ICE-99-6353	Coregel SemiPrep 87H	SS	21.2mm	300mm
<b>Glyphosate Analysis Columns:</b>				
GLY-99-2515	Glyphosate Guard Cartridge 2/pk	Pktit	4.0mm	20mm
GLY-99-3515	Glyphosate Guard Kit (Holder + 2 Cartridges)	SS/pktit	4.0mm	24mm
GLY-99-4610	Glyphosate Analysis Column, 100x4.6	SS	4.6mm	100mm
GLY-99-4630	Glyphosate Analysis Column, 300x4.6	SS	4.6mm	300mm
GLY-99-6515	Glyphosate Analysis Column, 150x4.6	SS	4.6mm	150mm
<b>Reversed Phase Columns:</b>				
AXC-99-1300	Universal Holder For All 4.0x20mm Guard Cartridges	SS	4.0mm	24mm
RPC-99-1314	RPSep PRX-1 Guard Cartridge 2/pk	Pktit	4.0mm	20mm
RPC-99-2314	RPSep PRX-1 Guard Kit (Holder + 2 Cartridges)	SS/pktit	4.0mm	24mm
RPC-99-7514	RPSep PRX-1 Column	SS	4.6mm	150mm
RPC-99-8514	RPSep PRX-1 Column	SS	4.6mm	250mm